

Product datasheet for **HP211759**

Adractalin (AAAS) Human qPCR Primer Pair (NM_015665)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	8086
Forward Sequence:	CCTTAGTGCCTCTGTCTTGGCT
Reverse Sequence:	GCTGGTAACAGGTGTATGCCCA
ACCN:	NM_015665 , NM_015665.1 , NM_015665.2 , NM_015665.3 , NM_015665.4 , NM_015665.5 , BC000659 , BC000659.2 , BC018836 , BT006912
UniProt ID:	Q9NRG9
Synonyms:	AAA; AAASb; ADRACALA; ADRACALIN; ALADIN; GL003
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)