

Product datasheet for **HP211812**

WT1-AS Human qPCR Primer Pair (NM_015855)

Product data:

| | |
|-------------------|--|
| Product Type: | qPCR Primer Pairs |
| Gene ID: | 51352 |
| Forward Sequence: | CCGGCCTCGAAGGTGCATAATT |
| Reverse Sequence: | CGCTCTGTCAGATGCAGTGGAC |
| ACCN: | BC002734 , NM_015855 , NM_015855.1 , NM_015855.2 , NM_015855.3 , BC096708 , BC098290 , BC100673 , BM543913 , BX110205 , BX387533 , BX387534 |
| Synonyms: | WIT-1, dj74j1.1, MGC120207, MGC120208, MGC120209 |
| Components: | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions) |
| Quality Control: | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage: | Store at -20°C. |
| Stability: | The primer mix is stable for one year from date of shipping. |



[View online »](#)