

## Product datasheet for **HP211836**

### EDD (UBR5) Human qPCR Primer Pair (NM\_015902)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	51366
Forward Sequence:	CGTGTCTACAGGACTGGAATGC
Reverse Sequence:	AACAGTCCAGCCGAATTGTGCC
ACCN:	<a href="#">NM_015902</a> , <a href="#">NM_015902.1</a> , <a href="#">NM_015902.2</a> , <a href="#">NM_015902.3</a> , <a href="#">NM_015902.4</a> , <a href="#">NM_015902.5</a> , <a href="#">BC037843</a> , <a href="#">BC044598</a> , <a href="#">BC053948</a> , <a href="#">BC137234</a> , <a href="#">BC151231</a> , <a href="#">BX489192</a> , <a href="#">NM_015902.6</a>
UniProt ID:	<a href="#">O95071</a>
Synonyms:	DD5; EDD; EDD1; HYD
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)