

Product datasheet for **HP211852**

Lipolysis Stimulated Lipoprotein Receptor (LSR) Human qPCR Primer Pair (NM_015925)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	51599
Forward Sequence:	GACGTTGACAGGAGTAGCTCAG
Reverse Sequence:	CCTTCTCCATGTAGTACAGGACC
ACCN:	BC004381 , NM_015925 , NM_015925.1 , NM_015925.2 , NM_015925.3 , NM_015925.5 , NM_015925.6 , BC004381.1 , BC047376 , BC047376.1 , BC000015 , BM837160 , BM971691 , BQ711340
UniProt ID:	Q86X29
Synonyms:	ILDR3; LISCH7
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)