

Product datasheet for **HP211853**

TEX264 Human qPCR Primer Pair (NM_015926)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	51368
Forward Sequence:	CATCCTGCCTTGGACACCTACA
Reverse Sequence:	CTCAGGCACATAGAAGTCTCCC
ACCN:	NM_015926 , NM_015926.1 , NM_015926.2 , NM_015926.3 , NM_015926.4 , NM_015926.5 , BC008742 , BC008742.2 , BC001592 , BG326586 , BG772477 , BI826172 , BP252919 , NM_015926.6
UniProt ID:	Q9Y6I9
Synonyms:	ZSIG11
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)