

Product datasheet for **HP211854**

LIPT1 Human qPCR Primer Pair (NM_015929)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	51601
Forward Sequence:	CAGGAACAGCTTCTAAGATCGGC
Reverse Sequence:	AGCAGTGGCATTGCTCCTGATC
ACCN:	NM_015929 , NM_015929.1 , NM_015929.2 , NM_015929.3 , BC070145 , BC070145.1 , BC007001 , BC009772 , BG717409 , BQ422048 , BQ432795 , BX100979 , NM_015929.4
UniProt ID:	Q9Y234
Synonyms:	LIPT1D
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)