

Product datasheet for **HP211871**

MRPL2 Human qPCR Primer Pair (NM_015950)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	51069
Forward Sequence:	GAGGAGAAGGTTATCCAAGTCCG
Reverse Sequence:	GATTGTATCTCCAGCCTGCATGT
ACCN:	BC020212 , NM_015950 , NM_015950.1 , NM_015950.2 , NM_015950.3 , NM_015950.4 , BC013685 , BM563754 , BM970997 , NM_015950.5
UniProt ID:	Q5T653
Synonyms:	CGI-22; MRP-L14; RPML14
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)