

Product datasheet for **HP211927**

MRPS2 Human qPCR Primer Pair (NM_016034)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	51116
Forward Sequence:	GTTTTCCGTGAGAAGCCTCTTCG
Reverse Sequence:	TGGCTGTCTGTTCCAGGTCGAT
ACCN:	NM_016034 , NM_016034.1 , NM_016034.2 , NM_016034.3 , NM_016034.4 , BC008017 , BC008017.1 , BC004905 , BC013108 , BE264885 , BI260471 , BQ065672 , BQ723448 , BU169388 , BX362434 , NM_016034.5
UniProt ID:	Q9Y399
Synonyms:	CGI-91; COXPD36; MRP-S2; S2mt
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)