

## Product datasheet for **HP211932**

### TMED5 Human qPCR Primer Pair (NM\_016040)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	50999
Forward Sequence:	AGAAGGAGTGCTTCTACCAGCC
Reverse Sequence:	CTAAGGTTTTGCCTTCTGGAGAG
ACCN:	<a href="#">NM_016040</a> , <a href="#">NM_016040.2</a> , <a href="#">NM_016040.3</a> , <a href="#">NM_016040.4</a> , <a href="#">BC070051</a> , <a href="#">BC070051.1</a> , <a href="#">BC016365</a> , <a href="#">BC016556</a> , <a href="#">BC023274</a> , <a href="#">BC038511</a> , <a href="#">BX642274</a>
UniProt ID:	<a href="#">Q9Y3A6</a>
Synonyms:	CGI-100; p24g2; p28
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)