

Product datasheet for **HP211949**

CIAO2B Human qPCR Primer Pair (NM_016062)

Product data:

| | |
|-------------------|--|
| Product Type: | qPCR Primer Pairs |
| Gene ID: | 51647 |
| Forward Sequence: | GATTCGCTCCATCAATGACCCG |
| Reverse Sequence: | CAGTGCGGAATGGTTGGTGTGA |
| ACCN: | NM_016062 , NM_016062.1 , NM_016062.2 , NM_016062.3 , BC005023 , BC005023.1 , BC001733 , BI598300 , BI603322 , BM091095 , BQ955804 , NM_016062.4 |
| UniProt ID: | Q9Y3D0 |
| Synonyms: | CGI-128; CIA2B; FAM96B; MIP18 |
| Components: | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions) |
| Quality Control: | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage: | Store at -20°C. |
| Stability: | The primer mix is stable for one year from date of shipping. |



[View online »](#)