

Product datasheet for **HP211950**

HDDC2 Human qPCR Primer Pair (NM_016063)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	51020
Forward Sequence:	TGTACGCCTAGCCCTGGTTCAT
Reverse Sequence:	TGGTAGGAGCTGGGTTATCTGC
ACCN:	BC003357 , NM_016063 , NM_016063.1 , BC003357.1 , BC066332 , BC066332.1 , BC001671 , BC061520 , BE743302 , BP257644
UniProt ID:	Q7Z4H3
Synonyms:	C6orf74; CGI-130; dJ167O5.2; NS5ATP2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)