

Product datasheet for **HP211971**

ZNF589 Human qPCR Primer Pair (NM_016089)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	51385
Forward Sequence:	TGGCTGTGCTTTTCACTGAGGC
Reverse Sequence:	AAGGGCAGGTATGGACTTCTGG
ACCN:	NM_016089 , NM_016089.1 , NM_016089.2 , BC005859 , BC006247 , BC008213 , BC028160 , BC048798 , BC062690 , BC166625 , NM_016089.3
UniProt ID:	Q86UQ0
Synonyms:	SZF1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)