

## Product datasheet for **HP211976**

### PSF2 (GINS2) Human qPCR Primer Pair (NM\_016095)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	51659
Forward Sequence:	AGCCAAACTCCGAGTGTCTGCT
Reverse Sequence:	CTTGTGTGAGGAAAGTCCCGCT
ACCN:	<a href="#">NM_016095</a> , <a href="#">NM_016095.1</a> , <a href="#">NM_016095.2</a> , <a href="#">BC003186</a> , <a href="#">BC003186.1</a> , <a href="#">BC010164</a> , <a href="#">BC022839</a> , <a href="#">BC062444</a> , <a href="#">BE514978</a> , <a href="#">NM_016095.3</a>
UniProt ID:	<a href="#">Q9Y248</a>
Synonyms:	HSPC037; Pfs2; PSF2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)