

Product datasheet for **HP211977**

IER3IP1 Human qPCR Primer Pair (NM_016097)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	51124
Forward Sequence:	CGAGGAGCGATTCTCAAGAAC
Reverse Sequence:	CTCATCACGGTTCTTACAGATCG
ACCN:	NM_016097 , NM_016097.1 , NM_016097.2 , NM_016097.3 , NM_016097.4 , BC010888 , BC010888.1 , BC017391 , BM829309 , NM_016097.5
UniProt ID:	Q9Y5U9
Synonyms:	HSPC039; MEDS; PRO2309
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)