

Product datasheet for **HP211978**

BRP44L (MPC1) Human qPCR Primer Pair (NM_016098)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	51660
Forward Sequence:	CCACGCAACAAATGAAGTAGCCC
Reverse Sequence:	CTGGCAATGCTGTCCCTCAAG
ACCN:	<u>NM_016098</u> , <u>NM_016098.1</u> , <u>NM_016098.2</u> , <u>BC000810</u> , <u>BC000810.1</u> , <u>BF665465</u> , <u>BG828174</u> , <u>NM_016098.4</u>
UniProt ID:	<u>Q9Y5U8</u>
Synonyms:	BRP44L; CGI-129; MPYCD; SLC54A1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)