

Product datasheet for **HP211981**

FKBP7 Human qPCR Primer Pair (NM_016105)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	51661
Forward Sequence:	AGAGGAGAGCACCGAAGAAGTG
Reverse Sequence:	TTTCGAGCCGTCTTTAGCCAGG
ACCN:	<u>NM_016105</u> , <u>NM_016105.2</u> , <u>BC009711</u> , <u>BC009711.1</u> , <u>BF210667</u> , <u>BQ003190</u> , <u>BT007122</u>
UniProt ID:	<u>Q9Y680</u>
Synonyms:	FK506 binding protein 7; FK506-binding protein 7; FK506-binding protein 23; FKBP23; FKBP23, PPIase, MGC9420; MGC9420; peptidyl-prolyl cis-trans isomerase; PPIase; rotamase
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)