

Product datasheet for **HP211993**

RHD Human qPCR Primer Pair (NM_016124)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	6007
Forward Sequence:	CAACGATACCCAGTTTGTCTGCC
Reverse Sequence:	GCATAGTAGGTGTTGAACACGGC
ACCN:	<u>NM_016124</u> , <u>NM_016124.1</u> , <u>NM_016124.2</u> , <u>NM_016124.3</u> , <u>NM_016124.4</u> , <u>BC139922</u> , <u>BC160075</u> , <u>BF798211</u>
UniProt ID:	<u>Q02161</u>
Synonyms:	CD240D; DIIIc; RH; Rh4; RH30; RHCED; RhDCw; RHDel; RHDVA(TT); RhII; RhK562-II; RhPI; RHPII; RHXIII
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)