

Product datasheet for **HP212007**

NSFL1C Human qPCR Primer Pair (NM_016143)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	55968
Forward Sequence:	CGACGAATCAGAGCCTACCACA
Reverse Sequence:	GCATCCACGATGAAGAGTCGGA
ACCN:	<u>NM_016143</u> , <u>NM_016143.2</u> , <u>NM_016143.3</u> , <u>NM_016143.4</u> , <u>BC002801</u> , <u>BC002801.2</u> , <u>BE795675</u> , <u>BQ007132</u>
UniProt ID:	<u>Q9UNZ2</u>
Synonyms:	dj776F14.1; P47; UBX1; UBXD10; UBXN2C
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)