

## Product datasheet for **HP212016**

### MMP17 Human qPCR Primer Pair (NM\_016155)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	4326
Forward Sequence:	CACCAAGTGAACAAGAGGAACC
Reverse Sequence:	GCAATGTCGCTCCAGACCTGA
ACCN:	<a href="#">NM_016155</a> , <a href="#">NM_016155.1</a> , <a href="#">NM_016155.2</a> , <a href="#">NM_016155.3</a> , <a href="#">NM_016155.4</a> , <a href="#">NM_016155.5</a> , <a href="#">BC045610</a> , <a href="#">BC045610.1</a> , <a href="#">BC040507</a> , <a href="#">BQ227546</a> , <a href="#">BX389435</a> , <a href="#">NM_016155.7</a>
UniProt ID:	<a href="#">Q9ULZ9</a>
Synonyms:	MMP-17; MT4-MMP; MT4MMP; MTMMP4
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)