

Product datasheet for **HP212017**

MTMR2 Human qPCR Primer Pair (NM_016156)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	8898
Forward Sequence:	ACGAGGAACTCTGACTGTCACG
Reverse Sequence:	TCACCTCGACTAGAAGCACCAC
ACCN:	NM_016156 , NM_016156.1 , NM_016156.2 , NM_016156.3 , NM_016156.4 , NM_016156.5 , BC052990 , BC052990.1 , BC012086 , BC034474 , BC040432 , BF960838 , BG252116 , BM475570 , BP358035 , BX107577 , BX647997
UniProt ID:	Q13614
Synonyms:	CMT4B; CMT4B1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)