

Product datasheet for **HP212029**

TRPC4 Human qPCR Primer Pair (NM_016179)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	7223
Forward Sequence:	CTTGAACAGGCAAGGTCCACCA
Reverse Sequence:	AAGTCCGCCATCCCACATCTGT
ACCN:	NM_016179 , NM_016179.1 , NM_016179.2 , BC104725 , BC105306 , BC105641 , BC105917 , BX648312 , NM_016179.4
UniProt ID:	Q9UBN4
Synonyms:	HTRP-4; HTRP4; TRP4
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)