

Product datasheet for **HP212031**

C1orf33 (MRTO4) Human qPCR Primer Pair (NM_016183)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	51154
Forward Sequence:	ACAGCAAGCTGAAGGACATCCG
Reverse Sequence:	GACCTGGTGCAGGTTGTCTTTG
ACCN:	<u>NM_016183</u> , <u>NM_016183.1</u> , <u>NM_016183.2</u> , <u>NM_016183.3</u> , <u>BC006504</u> , <u>BC006504.2</u> , <u>BC003013</u> , <u>BP312953</u> , <u>BQ670859</u>
UniProt ID:	<u>Q9UKD2</u>
Synonyms:	C1orf33; dj657E11.4; MRT4
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)