

Product datasheet for **HP212034**

BAF53b (ACTL6B) Human qPCR Primer Pair (NM_016188)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	51412
Forward Sequence:	CTCCACATCGACACCAATGCC
Reverse Sequence:	CAGGTTTGCTCAGACTTGACG
ACCN:	NM_016188 , NM_016188.1 , NM_016188.2 , NM_016188.3 , NM_016188.4 , BC020944 , BC020944.1 , BC008582 , BM696009 , NM_016188.5
UniProt ID:	O94805
Synonyms:	ACTL6; arpNalpha; BAF53B; DEE76; EIEE76; IDDSSAD
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)