

Product datasheet for **HP212037**

KIF20B Human qPCR Primer Pair (NM_016195)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	9585
Forward Sequence:	GCTGACTTTAAGGAGACTCTGCT
Reverse Sequence:	GTGGCACAATGTCTTTCGCTGC
ACCN:	NM_016195 , NM_016195.1 , NM_016195.2 , NM_016195.3 , BC005221 , BC012585 , BC042130 , BC046134 , BC058913 , BC093089 , BC108688 , NM_016195.4
UniProt ID:	Q96Q89
Synonyms:	CT90; KRMP1; MPHOSPH1; MPP-1; MPP1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)