

Product datasheet for **HP212043**

PDGFC Human qPCR Primer Pair (NM_016205)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	56034
Forward Sequence:	TGAACCAGGGTTCTGCATCCAC
Reverse Sequence:	TAAGCAGGTCCAGTGGCAAAGC
ACCN:	NM_016205 , NM_016205.1 , NM_016205.2 , BC041783 , BC051876 , BC094746 , BC136662 , BC144348 , NM_016205.3
UniProt ID:	Q9NRA1
Synonyms:	FALLOTEIN; SCDGF
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)