

## Product datasheet for **HP212087**

### LEF1 Human qPCR Primer Pair (NM\_016269)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	51176
Forward Sequence:	CTACCCATCCTCACTGTCTAGTC
Reverse Sequence:	GGATGTTCTGTTTGACCTGAGG
ACCN:	<a href="#">NM_016269</a> , <a href="#">NM_016269.1</a> , <a href="#">NM_016269.2</a> , <a href="#">NM_016269.3</a> , <a href="#">NM_016269.4</a> , <a href="#">BC050632</a> , <a href="#">BC050632.2</a> , <a href="#">BC040559</a> , <a href="#">NM_016269.5</a>
UniProt ID:	<a href="#">Q9UJU2</a>
Synonyms:	LEF-1; TCF1ALPHA; TCF7L3; TCF10
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)