

Product datasheet for **HP212097**

TAOK3 Human qPCR Primer Pair (NM_016281)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	51347
Forward Sequence:	CTACCCAGTCAGTTCAGAGCCA
Reverse Sequence:	CCTGACATCTGTTCCCGCAACT
ACCN:	NM_016281 , NM_016281.1 , NM_016281.2 , NM_016281.3 , BC002756 , BC002756.2 , BT007185 , BX649175 , NM_016281.4
UniProt ID:	Q9H2K8
Synonyms:	DPK; hKFC-A; JIK; MAP3K18
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)