

## Product datasheet for **HP212104**

### Hsp75 (TRAP1) Human qPCR Primer Pair (NM\_016292)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10131
Forward Sequence:	CTCTGTGGAGACGGACATAGTC
Reverse Sequence:	CTCCTTCTCTGATAGGCACTCG
ACCN:	<a href="#">NM_016292</a> , <a href="#">NM_016292.1</a> , <a href="#">NM_016292.2</a> , <a href="#">BC018950</a> , <a href="#">BC000406</a> , <a href="#">BC001455</a> , <a href="#">BC002994</a> , <a href="#">BC023585</a> , <a href="#">NM_016292.3</a>
UniProt ID:	<a href="#">Q12931</a>
Synonyms:	HSP 75; HSP75; HSP90L; TRAP-1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)