

Product datasheet for **HP212129**

EVL Human qPCR Primer Pair (NM_016337)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	51466
Forward Sequence:	GGACCTCAAAGTCCGATGCCAA
Reverse Sequence:	TGTCTGACTGGGAGGCTGCTTT
ACCN:	NM_016337 , NM_016337.1 , NM_016337.2 , BC032358 , BC032358.1 , BC023997 , BC049376 , BQ083568 , NM_016337.3
UniProt ID:	Q9UI08
Synonyms:	RNB6
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)