

## Product datasheet for **HP212134**

### CENPF Human qPCR Primer Pair (NM\_016343)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	1063
Forward Sequence:	AGCACGACTCCAGCTACAAGGT
Reverse Sequence:	CATCATGCTTTGGTGTCTTTCTG
ACCN:	<a href="#">NM_016343</a> , <a href="#">NM_016343.1</a> , <a href="#">NM_016343.2</a> , <a href="#">NM_016343.3</a> , <a href="#">BC172232</a> , <a href="#">BG680819</a> , <a href="#">BQ774943</a> , <a href="#">BU155312</a> , <a href="#">BU662228</a> , <a href="#">NM_016343.4</a>
UniProt ID:	<a href="#">P49454</a>
Synonyms:	CENF; CILD31; hcp-1; PRO1779; STROMS
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)