

Product datasheet for **HP212143**

DCDC2 Human qPCR Primer Pair (NM_016356)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	51473
Forward Sequence:	TGACAAGTCAACGATGAGAAGGC
Reverse Sequence:	CCAACTGTTGACTTAGAGTGGCG
ACCN:	<u>NM_016356</u> , <u>NM_016356.1</u> , <u>NM_016356.2</u> , <u>NM_016356.3</u> , <u>NM_016356.4</u> , <u>BC050704</u> , <u>BC050704.1</u> , <u>BF109922</u> , <u>BM923267</u> , <u>NM_016356.5</u>
UniProt ID:	<u>Q9UHG0</u>
Synonyms:	DCDC2A; DFNB66; NPHP19; NSC; RU2; RU2S
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)