

Product datasheet for **HP212153**

TPRA1 Human qPCR Primer Pair (NM_016372)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	131601
Forward Sequence:	GTGTAGATGCCACAACCTTCCTG
Reverse Sequence:	CTCTGTCTCGTCCACTTGGCAT
ACCN:	<u>NM_016372</u> , <u>NM_016372.1</u> , <u>NM_016372.2</u> , <u>BC050711</u> , <u>BC050711.2</u> , <u>BC017058</u> , <u>BG168441</u> , <u>BI823558</u> , <u>BM665705</u>
UniProt ID:	<u>Q86W33</u>
Synonyms:	GPR175; TMEM227; TPRA40
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)