

Product datasheet for **HP212154**

WVOX Human qPCR Primer Pair (NM_016373)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	51741
Forward Sequence:	TCCTCCAGGATGTTTTGTGCCG
Reverse Sequence:	AAGCCAGCATCGCCAATAGTC
ACCN:	<u>NM_016373</u> , <u>NM_016373.1</u> , <u>NM_016373.2</u> , <u>NM_016373.3</u> , <u>BC003184</u> , <u>BT007445</u> , <u>BU157392</u> , <u>NM_016373.4</u>
UniProt ID:	<u>Q9NZC7</u>
Synonyms:	D16S432E; DEE28; EIEE28; FOR; FRA16D; HHCMA56; PRO0128; SCAR12; SDR41C1; WOX1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)