

## Product datasheet for **HP212179**

### RNF141 Human qPCR Primer Pair (NM\_016422)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	50862
Forward Sequence:	TCCTCTGTAACATCTTGTCAGGC
Reverse Sequence:	CTGTGAGCACAAAGGCAGGATGA
ACCN:	<a href="#">NM_016422</a> , <a href="#">NM_016422.1</a> , <a href="#">NM_016422.2</a> , <a href="#">NM_016422.3</a> , <a href="#">BC018104</a> , <a href="#">BC018104.1</a> , <a href="#">BC035089</a> , <a href="#">BF970569</a> , <a href="#">BG108453</a> , <a href="#">BQ229248</a> , <a href="#">BT006662</a> , <a href="#">BU679024</a> , <a href="#">BX647862</a>
UniProt ID:	<a href="#">Q8WVD5</a>
Synonyms:	RFP141; ZFP26; ZNF230
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)