

## Product datasheet for **HP212180**

### GTSE1 Human qPCR Primer Pair (NM\_016426)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	51512
Forward Sequence:	CTCTACCAGCAATCTCGCAAGG
Reverse Sequence:	GACTTGCTGATGTTTGACAGAGG
ACCN:	<u><a href="#">BC021234</a></u> , <u><a href="#">NM_016426</a></u> , <u><a href="#">NM_016426.1</a></u> , <u><a href="#">NM_016426.2</a></u> , <u><a href="#">NM_016426.3</a></u> , <u><a href="#">NM_016426.4</a></u> , <u><a href="#">NM_016426.5</a></u> , <u><a href="#">NM_016426.6</a></u> , <u><a href="#">BC006325</a></u> , <u><a href="#">BC069212</a></u> , <u><a href="#">BC075828</a></u> , <u><a href="#">BE390871</a></u> , <u><a href="#">BG718574</a></u> , <u><a href="#">BG770679</a></u> , <u><a href="#">BQ049956</a></u> , <u><a href="#">BU621031</a></u> , <u><a href="#">BX329502</a></u>
UniProt ID:	<u><a href="#">Q9NYZ3</a></u>
Synonyms:	B99
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)