

Product datasheet for **HP212189**

HIGD1B Human qPCR Primer Pair (NM_016438)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	51751
Forward Sequence:	TGAGGTCTCGTGGTTCCACCAA
Reverse Sequence:	GCCATCCTCTTGACGTAATCGC
ACCN:	NM_016438 , NM_016438.1 , NM_016438.2 , NM_016438.3 , BC020667 , BC020667.1 , BE041447 , BG573492 , BM503772 , BM876435 , BQ017345 , BQ292861 , NM_016438.4
UniProt ID:	Q9P298
Synonyms:	CLST11240; CLST11240-15
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)