

Product datasheet for **HP212193**

CDT2 (DTL) Human qPCR Primer Pair (NM_016448)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	51514
Forward Sequence:	CCAGCCTTAGTCCAGATGACCA
Reverse Sequence:	GAGAATGACCCAGGAGCACAGT
ACCN:	<u>NM_016448</u> , <u>NM_016448.1</u> , <u>NM_016448.2</u> , <u>NM_016448.3</u> , <u>BC033540</u> , <u>BC033540.1</u> , <u>BC033297</u> , <u>BF508267</u> , <u>NM_016448.4</u>
UniProt ID:	<u>Q9NZJ0</u>
Synonyms:	CDT2; DCAF2; L2DTL; RAMP
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)