

Product datasheet for **HP212195**

beta COP (COPB1) Human qPCR Primer Pair (NM_016451)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	1315
Forward Sequence:	GCAACTCAGAGTGCCTTAGCA
Reverse Sequence:	GCAATCTTGGTCAGAGTTGTGGC
ACCN:	NM_016451 , NM_016451.1 , NM_016451.2 , NM_016451.3 , NM_016451.4 , BC037280 , BC037280.2 , BC012572
UniProt ID:	P53618
Synonyms:	BARMACS; COPB
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)