

## Product datasheet for **HP212229**

### MRPL27 Human qPCR Primer Pair (NM\_016504)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	51264
Forward Sequence:	CTCTTGCTGTCAGATACGCATCC
Reverse Sequence:	CGGAAATGGCGCTGTGTTGCAA
ACCN:	<a href="#">NM_016504</a> , <a href="#">NM_016504.1</a> , <a href="#">NM_016504.2</a> , <a href="#">BC001066</a> , <a href="#">BC001066.2</a> , <a href="#">BC021986</a> , <a href="#">BC021986.1</a> , <a href="#">BG829493</a> , <a href="#">BM840463</a> , <a href="#">NM_016504.3</a>
UniProt ID:	<a href="#">Q9P0M9</a>
Synonyms:	L27mt
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)