

## Product datasheet for **HP212232**

### SCLY Human qPCR Primer Pair (NM\_016510)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	51540
Forward Sequence:	ACTTCCTCGGTGGAACACGACT
Reverse Sequence:	ATGATGGTCACGAGCGTGTGG
ACCN:	<u><a href="#">NM_016510</a></u> , <u><a href="#">NM_016510.1</a></u> , <u><a href="#">NM_016510.2</a></u> , <u><a href="#">NM_016510.3</a></u> , <u><a href="#">NM_016510.4</a></u> , <u><a href="#">NM_016510.5</a></u> , <u><a href="#">BC007891</a></u> , <u><a href="#">BC007891.2</a></u> , <u><a href="#">BC000586</a></u> , <u><a href="#">BC018966</a></u>
UniProt ID:	<u><a href="#">Q96I15</a></u>
Synonyms:	hSCL; SCL
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)