

Product datasheet for **HP212250**

SIRT6 Human qPCR Primer Pair (NM_016539)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	51548
Forward Sequence:	TGGCAGTCTTCCAGTGTGGTGT
Reverse Sequence:	CGCTCTCAAAGGTGGTGTGAA
ACCN:	<u>BC005026</u> , <u>NM_016539</u> , <u>NM_016539.1</u> , <u>NM_016539.2</u> , <u>NM_016539.3</u> , <u>BC005026.1</u> , <u>BC028220</u> , <u>BC004218</u> , <u>BE281103</u> , <u>BP418952</u> , <u>BQ060703</u> , <u>NM_016539.4</u>
UniProt ID:	<u>Q8N6T7</u>
Synonyms:	SIR2L6
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)