

Product datasheet for **HP212263**

BFAR Human qPCR Primer Pair (NM_016561)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	51283
Forward Sequence:	AGGCTGAGTCTGCTCTACCTGT
Reverse Sequence:	ACTGCTTCCACGTAGGCTCCTT
ACCN:	<u>NM_016561</u> , <u>NM_016561.1</u> , <u>NM_016561.2</u> , <u>BC003054</u> , <u>BC003054.1</u> , <u>BI912683</u> , <u>BU621058</u> , <u>NM_016561.3</u>
UniProt ID:	<u>Q9NZS9</u>
Synonyms:	BAR; RNF47
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)