

Product datasheet for **HP212271**

LGSN Human qPCR Primer Pair (NM_016571)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	51557
Forward Sequence:	GGTCTAAGACTATCCCTGCACAC
Reverse Sequence:	CTCTGGCATTAGGACTATGTCGC
ACCN:	<u>NM_016571</u> , <u>NM_016571.1</u> , <u>NM_016571.2</u> , <u>BC130366</u> , <u>BC130368</u> , <u>NM_016571.3</u>
UniProt ID:	<u>Q5TDP6</u>
Synonyms:	GLULD1; LGS
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)