

Product datasheet for **HP212324**

CRNKL1 Human qPCR Primer Pair (NM_016652)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	51340
Forward Sequence:	CCATTCAGGAGAAGAGGCACTG
Reverse Sequence:	GGCTTGATACACCTGTCTTGTC
ACCN:	NM_016652 , NM_016652.1 , NM_016652.2 , NM_016652.3 , NM_016652.4 , NM_016652.5 , BC144222 , BC156789 , BI550198
UniProt ID:	Q9BZJ0
Synonyms:	CLF; Clf1; CRN; HCRN; MSTP021; SYF3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)