

Product datasheet for **HP212328**

SHC3 Human qPCR Primer Pair (NM_016848)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	53358
Forward Sequence:	CACTTTTGCGAAGACTGGCAG
Reverse Sequence:	CTGCTGAGTGTGACGTAGGTG
ACCN:	NM_016848 , NM_016848.1 , NM_016848.2 , NM_016848.3 , NM_016848.4 , NM_016848.5 , BC026314 , BC026314.1 , BC042431 , BX538228 , BX641139
UniProt ID:	Q92529
Synonyms:	N-Shc; NSHC; RAI; SHCC
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)