

Product datasheet for **HP212331**

SIX2 Human qPCR Primer Pair (NM_016932)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10736
Forward Sequence:	CACACAGGTCAGCAACTGGTTC
Reverse Sequence:	TCATCCTCCGAGCTGCCTAACA
ACCN:	NM_016932 , NM_016932.1 , NM_016932.2 , NM_016932.3 , NM_016932.4 , BC024033 , BC024033.2 , BT020100
UniProt ID:	Q9NPC8
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)