

Product datasheet for **HP212342**

Homeo box C10 (HOXC10) Human qPCR Primer Pair (NM_017409)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	3226
Forward Sequence:	GAGCGAAAAGGAGAGGGCCAAA
Reverse Sequence:	TCCGCTCTTTGCTGTCAGCCAA
ACCN:	<u>NM_017409</u> , <u>NM_017409.1</u> , <u>NM_017409.2</u> , <u>NM_017409.3</u> , <u>BC001293</u> , <u>BC001293.1</u> , <u>BT009921</u> , <u>NM_017409.4</u>
UniProt ID:	<u>Q9NYD6</u>
Synonyms:	HOX3I
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)