

Product datasheet for **HP212347**

UBP43 (USP18) Human qPCR Primer Pair (NM_017414)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	11274
Forward Sequence:	TGGACAGACCTGCTGCCTTAAC
Reverse Sequence:	CTGTCCTGCATCTTCTCCAGCA
ACCN:	<u>NM_017414</u> , <u>NM_017414.1</u> , <u>NM_017414.2</u> , <u>NM_017414.3</u> , <u>BC014896</u> , <u>BC014896.1</u> , <u>BI752849</u> , <u>BT006835</u> , <u>BU584005</u> , <u>NM_017414.4</u>
UniProt ID:	<u>Q9UMW8</u>
Synonyms:	ISG43; PTORCH2; UBP43
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)