

Product datasheet for **HP212366**

A4GALT Human qPCR Primer Pair (NM_017436)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	53947
Forward Sequence:	GCATCTACCTGGACACGGACTT
Reverse Sequence:	ATGCACAGCGCCATGAACTCGT
ACCN:	<u>BC017068</u> , <u>NM_017436</u> , <u>NM_017436.1</u> , <u>NM_017436.2</u> , <u>NM_017436.3</u> , <u>NM_017436.4</u> , <u>NM_017436.5</u> , <u>NM_017436.6</u> , <u>BC017068.2</u> , <u>BC055286</u> , <u>NM_017436.7</u>
UniProt ID:	<u>Q9NPC4</u>
Synonyms:	A4GALT1; A14GALT; Gb3S; P(k); P1; P1PK; PK
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)