

Product datasheet for **HP212382**

Kallikrein 15 (KLK15) Human qPCR Primer Pair (NM_017509)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	55554
Forward Sequence:	TAACGTGGCGCTTCCCTCATC
Reverse Sequence:	AGACGTGGTCCGTAGTTGCTCT
ACCN:	<u>NM_017509</u> , <u>NM_017509.1</u> , <u>NM_017509.2</u> , <u>NM_017509.3</u> , <u>BC069507</u> , <u>BC069507.1</u> , <u>BC069480</u> , <u>BC069518</u> , <u>BC126137</u> , <u>BC144045</u> , <u>BC144046</u> , <u>BC144047</u>
UniProt ID:	<u>Q9H2R5</u>
Synonyms:	ACO; HSRNASPH
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)